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## Which cola is the strongest ?

The aim of this work is to discover informations about the pH, its measurements, then to measure the pH of some solutions and to compare the pH of two sodas<sup>2</sup>.

## Documents : about pH measurements


## ACID OR BASIC ?

In chemistry, pH is a measure of the acidity or basicity of an aqueous solution. Pure water is said to be neutral, with a pH close to 7.  
Solutions with a pH inferior to 7 are said to be acidic and solutions with a pH superior to 7 are basic (or alkaline).  
pH measurements are important in medicine, biology, chemistry, agriculture, food science, environmental science, oceanography...

## PH SCALE

pH values are in a range of values between 0 to 14, from acid solutions (pH<7) to alkaline solutions (pH>7).

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

## USING PH STRIPS TO MEASURE THE PH OF A SOLUTION

|                                                                                   |                                                                                   |                                                                                    |                                                                                     |                                                                                     |
|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
|  |  |  |  |  |
| pH strips                                                                         | Step 1                                                                            | Step 2                                                                             | Step 3                                                                              | Step 4                                                                              |

## USING A PH METER TO MEASURE THE PH OF A SOLUTION

|                                                                                   |                                                                                   |                                                                                    |                                                                                     |                                                                                     |
|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
|  |  |  |  |  |
| pH meter                                                                          | Step 1                                                                            | Step 2                                                                             | Step 3                                                                              | Step 4                                                                              |

## Part one : Training

Measure the pH of the following solutions

|                                       | USING PH STRIPS | USING PH METER | ACIDIC, NEUTRAL OR BASIC ? |
|---------------------------------------|-----------------|----------------|----------------------------|
| <b>Solution 1 : hydrochloric acid</b> | 1               | 1,1            |                            |
| <b>Solution 2 : lemon juice</b>       |                 |                |                            |
| <b>Solution 3 : white vinegar</b>     |                 |                |                            |
| <b>Solution 4 : grapefruit juice</b>  |                 |                |                            |
| <b>Solution 5 : water</b>             |                 |                | neutral                    |
| <b>Solution 6 : milk</b>              |                 |                |                            |
| <b>Solution 7 : drain cleaner</b>     |                 |                |                            |
| <b>Solution 8 : soap</b>              |                 |                |                            |
| <b>Solution 9 : whitewash</b>         |                 |                |                            |

## Part two : which cola is the strongest ?

## USE THE DIFFERENT DOCUMENTS TO :

- discover which soda is the "strongest";
- discover if carbon dioxide is responsible of this acidity.

Your report should include your experience schemes, your observations and a conclusion.

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